

# MADomizer

## Bottle Atomizer



## The simple, effective and targeted topical medication delivery system\*

### Simple, Effective and Targeted

The MADomizer Bottle Atomizer from Teleflex is safe, effective and customizable to the individual needs of your patients.

The MADomizer Device allows exact dosing and targeted delivery of topical anesthetics and vasoconstrictors across a broad mucosal surface of the paranasal sinuses in a matter of seconds.

### Reduces the Risk of Cross-Contamination


The disposable tip and unique positive displacement pump reduce the risk of

cross-contamination typically associated with Venturi atomizers.


### Reliable, Flexible and Portable

The MADomizer Device is a reliable solution that has no parts to clean, does not clog, and allows for approximately 200 pumps per full bottle. The device is simple and easy to use, and has a malleable disposable cannula, which is not made with natural rubber latex, that allows for use with patients in any position. The device is portable, fitting on ENT carts and in consult bags, is individually packaged, and does not need noisy compressed air.

### Venturi Atomizer



### MADomizer with positive displacement pump



### Benefits



**The Clinician**  
Allows for exact dosing, is simple and easy to use



**Your Institution**  
Simple, effective and targeted topical medication delivery system that reduces the risk of cross contamination




**The Patient**  
Safe, effective and customizable to the individual needs of each patient

\*For use with drugs approved for intranasal and oropharyngeal delivery.



### Using the MADomizer Bottle Atomizer



- 1 Place the disposable applicator (MAD510) tip in nostril or oropharyngeal cavity
- 2 Compress the pump to spray solution into nostril or oropharyngeal cavity
- 3 Discard applicator (MAD510) after each patient use
- 4 Clean the exterior surface of the bottle and pump. Always select a hospital grade disinfectant wipe and use in accordance with the manufacturer's instructions
- 5 Place protective cap on pump stem for storage

Repeat steps 1–5 until all of the solution has been consumed or is no longer needed. The glass bottle and pump may not be reused. Dispose of glass bottle and pump in accordance with federal, state and local regulations.

#### MADomizer Device Specifications

Typical Particle Size	30–100 microns
Spray Volume	0.1 ml/spray
Tip Diameter	0.19" (4.82mm)
Applicator Length	4.75" (12.10 cm) end of tip to end of thumb landing / 2.25" (5.70 cm) end of tip to end of clip
Bottle Volume	20ml

#### MADomizer Bottle Atomizer Device Ordering

ITEM NUMBER	PRODUCT DESCRIPTION	QTY/BOX
MAD500	MADomizer Bottle Atomizer Device	5
MAD510	MADomizer Disposable Tip	50
MAD515	MADomizer Pull Tab Cap	20
MAD520	MADomizer Replacement Bottle Clip	5

To learn more, visit [Imaco.com/products](http://Imaco.com/products)

Teleflex, the Teleflex Logo, MAD, and MADomizer are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this document may not be available in all countries. Please contact your local representative. All data current at time of printing (05/2016). Subject to technical changes without further notice. © 2016 Teleflex Incorporated. All rights reserved. 94 07 99 - 00 00 01 - REV A - 05 15 01 - MC-000358 Rev. 1



Distributed by:  
**Teleflex Headquarters International, Ireland** · Teleflex Medical Europe Ltd. · IDA Business & Technology Park  
 Dublin Road · Athlone · Co Westmeath · Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · orders.intl@teleflex.com  
**United Kingdom** Tel. +44 (0)14 94 53 27 61 · info.uk@teleflex.com  
**South Africa** Tel. +27 (0)11 807 4887 · assist.africa@teleflex.com