


LMA MADdy™ Paediatric malleable topical anaesthetic device.



LMA MADdy™ is a single-use atomization device. Its unique design has a flexible cannula with memory and a tip that atomizes drugs into a fine mist of particles 30-100 microns in size.¹

Topical anaesthetic delivered to the mucosal surface of the larynx may encourage greater haemodynamic stability and reduced airway stimulation in patients:^{2,3}

- ▶ Reduced heartrate vs pre-induction values²
- ▶ Reduced cardiovascular response³
- ▶ Easier intubation vs saline²
- ▶ Fewer attempts for successful insertion of a laryngeal mask⁴

Laryngospasms are three times more likely to occur in paediatric patients than an adult.

27.6 Incidence of laryngospasm per 1000 paediatric patients⁶

8.6 Incidence of laryngospasm per 1000 adult patients⁶

“An ideal topical anesthetic should have rapid onset, adequate duration, and minimal or predictable systemic absorption, with the ability to be accurately administered.”

VHA Pharmacy Benefits Management Strategic Healthcare Group and the Medical Advisory Panel and the National Center for Patient Safety, 2006





Atomization spray

The spray atomizes drugs into a fine mist of particles 30-100 microns in size.⁷

Syringe

Available with or without a syringe attached.

Pressure

High applied pressure ensures that drugs are atomized into a fine mist of particles through the tip of the device.

Blowfish


The Blowfish acts as a distracting agent to improve patient compliance.

Malleable stylet

360° flexibility of the stylet enables it to be adapted to an individual patient's anatomy,⁸ ensuring the targeted delivery of medication.

Accurate dosing

The syringe enables the accurate measurement of drug to be delivered.




LMA MADdy™: Specifications

	MAD900	MAD800
Typical particle size	30-100 microns	
Tip diameter	4.82 mm (0.19 inches)	
Tip length	11.4 cm (4.5 inches)	11.4 cm (4.5 inches)
System dead space	0.13 ml	0.19 ml

LMA MADdy™: Product specification

Product code	Quantity/case	Description
MAD900	25	Paediatric sized mucosal atomization device for atomizing topical solutions without syringe.
MAD800	25	Paediatric sized mucosal atomization device for atomizing topical solutions with 3 ml syringe.

Distributed by:



30-0336 Rev 02/14

TELEFLEX HEADQUARTERS INTERNATIONAL, IRELAND
Teleflex Medical Europe Ltd., IDA Business and Technology Park,
Dublin Road, Athlone, Co Westmeath, Ireland
Phone +353 (0)9 06 46 08 00 Fax +353 (0)14 37 07 73
orders.intl@teleflex.com



Manufactured by:
Wolfe Tory Medical, Inc.
79 West 4500 South
Salt Lake City, Utah 84107



Consult IFU on this website:
www.LMACO.com



Printed on paper from sustainable sources

Teleflex www.teleflex.com LMA www.lmaco.com

Copyright © 2014 Teleflex Incorporated. All rights reserved. LMA, LMA MAD, LMA MADdy, MAD, MADdy and LMA Better by Design are trademarks or registered trademarks of Teleflex Incorporated or its affiliates. Wolfe-Tory Medical, Inc. is a wholly owned subsidiary of Teleflex Medical.

XX XX XX – XX XX XX . XXX XXXX . XX/X . XX XX XX

References:

- Joshi N., Marjot R. Anaesthesia 2008; **63**: 1255-1267.
- Bülow K. et al. Acta Anaesthesiol Scand 1996; **40**: 752-756.
- Hamaya Y., Dohi S. Anesthesiologist 2000; **93**: 95-103.
- Cook T.M. et al. Anaesthesia 1996; **51**: 787-790.
- VHA Pharmacy Benefits Management Strategic Healthcare Group and the Medical Advisory Panel and the National Center for Patient Safety, 2006.
- Olsson G.L., Hallen B. Acta Anaesthesiol Scand 1984; **28**: 567-575.
- Talon M. et al. J Burn Care Res 2009; **30**: 599-605.
- Xue F.S. et al. BJA 2007; **99** (6): 920-921.



Clinical evidence



For the latest clinical evidence on LMA MADdy™

www.lmaco.com/evidence



For the latest educational videos and clinician testimonials

www.youtube.com/LaryngealMaskAirway



For the latest news from LMA

www.facebook.com/LMAInternational



Better by Design



For product information and access to product instructions for use

www.lmaco.com



For information on other products within the Teleflex product portfolio

www.teleflex.com