



LMA Unique™

The original single use laryngeal mask.

LMA Unique™ is a first generation, versatile, single use airway with aperture bars used in millions of surgical cases each year.

- ▶ Seven sizes, from neonate to adult
- ▶ Not made with natural rubber latex, sterile, single use
- ▶ Aperture bars designed to prevent the blockage of the airway by the epiglottis¹
- ▶ Soft, flexible cuff
- ▶ Facilitates smooth emergence from anaesthesia²
- ▶ Minimal haemodynamic response²
- ▶ Part of the ASA difficult airway algorithm³

Consider using LMA Unique™ for:

- ▶ Theatre
- ▶ Crash carts
- ▶ Emergency departments
- ▶ Pre-hospital
- ▶ Ambulance



“Advantages over the TT included: increased speed and ease of placement by inexperienced personnel; increased speed of placement by anaesthetists; improved haemodynamic stability at induction and during emergence; minimal increase in intraocular pressure following insertion; reduced anaesthetic requirements for airway tolerance; lower frequency of coughing during emergence; improved oxygen saturation during emergence; and lower incidence of sore throat in adults.”

Brimacombe J., 1995




LMA Unique™: Product specification

Mask size	Product code	Patient size	Maximum cuff volume (air)*
1	125010	Neonates/infants up to 5 kg	4 ml
1.5	125015	Infants 5-10 kg	7 ml
2	125020	Infants/children 10-20 kg	10 ml
2.5	125025	Children 20-30 kg	14 ml
3	125030	Children 30-50 kg	20 ml
4	125040	Adults 50-70 kg	30 ml
5	125050	Adults 70-100 kg	40 ml

*These are maximum volumes that should never be exceeded. It is recommended that the cuff be inflated to a maximum of 60 cm H₂O intracuff pressure.

Distributed by:



PAF-2500-005 Rev B IW 201402

TELEFLEX HEADQUARTERS INTERNATIONAL, IRELAND
 Teleflex Medical Europe Ltd., IDA Business and Technology Park,
 Dublin Road, Athlone, Co Westmeath, Ireland
 Phone +353 (0)9 06 46 08 00 Fax +353 (0)14 37 07 73
 orders.intl@teleflex.com



Manufactured by:
 The Laryngeal Mask Company Limited
 Le Rocher, Victoria, Mahé, Seychelles



Consult IFU on this website:
www.LMACO.com



Printed on paper from sustainable sources

Teleflex www.teleflex.com LMA www.lmaco.com

Copyright © 2014 Teleflex Incorporated. All rights reserved. LMA, LMA Unique, LMA Better by Design are trademarks or registered trademarks of Teleflex Incorporated or its affiliates.

XX XX XX – XX XX XX . XXX XX/XX . XX/X . XX XX XX

References:

1. Benumof J.L. CJA 2003; **50**: 968.
2. Brimacombe J. CJA 1995; **42**: 1017-1023.
3. American Society of Anaesthesiologists Task Force on Management of the Difficult Airway. Anesthesiol 2003; **98**: 1269-1277.



Clinical evidence



For the latest clinical evidence on LMA Unique™

www.lmaco.com/evidence



YouTube



For the latest educational videos and clinician testimonials

www.youtube.com/LaryngealMaskAirway



For the latest news from LMA

www.facebook.com/LMAInternational



Better by Design



For product information and access to product instructions for use

www.lmaco.com



For information on other products within the Teleflex product portfolio

www.teleflex.com