

Silatos™ Silicone Sheeting



The Silatos™ Silicone Sheatings are made of a medical grade, flexible and translucent silicone elastomer material.

*Let's
talk about
life®*

ATOS
MEDICAL

The Silatos™ Silicone Sheetings are made of a medical grade, flexible and translucent silicone elastomer material.

The sheetings are easily cut to any size and shape. They can be autoclaved and are also highly permeable to gases. Removal of the sheetings is facilitated by their low tendency to adhere to adjacent tissue.

The Silatos™ Silicone Sheetings are designed to fit a variety of medical applications and are available in several thicknesses, flexibilities and with a double knit polyester fabric which reinforces the sheeting for dimensional stability and tear resistance.

Common Application

- Dressing of the ear canal after all kinds of middle ear surgery.
- Different applications in rhinoplastic surgery.
- Skin covering and treatment of cicatrices.
- Various laboratory uses.

Thorough testing is recommended for products made in part or otherwise incorporating Silatos Silicone Sheetings to determine the acceptability of the product's performance in a specific application.

For single use and for single patient use only!

PACKAGING

The sheetings are non-sterile and individually packed. The following designs are available:

| DESCRIPTION | REF | THICKNESS | SIZE |
|----------------|------|----------------------|--------------|
| Non-reinforced | 7458 | 0.13 mm (0.005 inch) | 150 x 200 mm |
| Non-reinforced | 7459 | 0.50 mm (0.020 inch) | 150 x 200 mm |
| Non-reinforced | 7457 | 1.00 mm (0.040 inch) | 150 x 200 mm |
| Reinforced | 7456 | 1.00 mm (0.040 inch) | 150 x 200 mm |
| Reinforced | 7455 | 0.50 mm (0.020 inch) | 150 x 200 mm |

Head office:
Atos Medical AB
 P.O. Box 183, SE-242 22
 Hörby, Sweden
 Tel: +46 (0)415 198 00
 Fax: +46 (0)415 198 98
 info@atosmedical.com
 www.atosmedical.com

U.K. Sales office:
 Platon Medical Ltd
 P O Box 2568, Eastbourne,
 East Sussex BN21 3HZ,
 United Kingdom
 Tel: 01323 431 930
 Fax: 01323 732 606
 platon@mistral.co.uk

Niederlassung:
 Atos Medical GmbH
 Carl Schuricht Strasse 3,
 DE-65187, Wiesbaden,
 Germany
 Tel: 0611 445 70-0
 Fax: 0611 445 7023
 info.de@atosmedical.com

Verkoopadres in Nederland:
 Atos Medical BV
 Blauw-roodlaan 286
 Postbus 574,
 2700 AN Zoetermeer
 The Netherlands
 Tel: 079 593 5000
 Fax: 079 593 5009
 Info.nl@atosmedical.com

Oficina de ventas en España:
 Atos Medical Spain S.L.
 C/Aragon 208-210 3º 3ª,
 ES-08011, Barcelona
 Spain
 Tel: 93 323 9195
 Fax: 93 451 3695
 info.es@atosmedical.com

**Bureau de vente en Belgique/
 Verkoopkantoor in België**
 Atos Medical BVBA / SPRL
 Koningsstraat 266 Rue Royale
 BE - 1210 Brussel / Bruxelles
 Tel 02 / 218 55 50
 Fax 02 / 218 55 51
 info.be@atosmedical.com

Atos Medical Switzerland
 Mainaustrasse 15
 Postfach
 8034 Zürich, Switzerland
 Tel: +41 (0)840 480 490
 Fax: +41 (0)840 480 400
 info.ch@atosmedical.com



© Atos Medical AB, 200610, REF 7940.