


EXTRA DEPENDABLE PRODUCTS FROM KAPITEX HEALTHCARE

# dysphagia management



**Featuring:**

- Kapi-Cup®
- Provale™ Cup
- Novo Cup
- Sure Grip Mug
- Dysphagia Cup
- Maroon Spoons



**KAPITEX**

The Specialists in Laryngectomy,  
Tracheostomy, Dysphagia and  
Airway Management

# KAPI-CUP®

Available in 3 sizes in packs of 5




These semi-flexible lightweight cups are carefully shaped to encourage safe drinking and swallowing - as the user does not need to tilt the head whilst drinking

- ✓ Semi-flexible top part of the cup allows the cup to be squeezed to help direct flow of fluid into the mouth by a carer or user
- ✓ Bite-safe and shatterproof - smooth moulded edges and lightweight
- ✓ Suitable for warm and cold liquids
- ✓ Competitively and economically priced
- ✓ Small size (30ml) ideal for use with children
- ✓ Each cup has an approximate fill line
- ✓ The medium and large sized cups are supplied with a double handle holder - making the price even more attractive!
- ✓ A choice of single or double handle holders are available in packs of 2
- ✓ Microwave and dishwasher safe
- ✓ CE marked under the Medical Device Directive 93/42/EEC



Single handle cup holder




Double handle cup holder

OS FSA 5501	Kapi-Cup Small	30ML	Green	Pack of 5
OS FSA 5502	Kapi-Cup Medium	80ML	Blue	Pack of 5*
OS FSA 5503	Kapi-Cup Large	160ML	White	Pack of 5*
OS FSA 5504	Kapi-Cup Holder Single Handle			Pack of 2
OS FSA 5505	Kapi-Cup Holder Double Handle			Pack of 2

\*Supplied with one double handle holder

# SURE GRIP MUG

Ergonomically designed, easy to hold, even in small or weak hands. With Sure Grip mugs the user has complete control over the flow of liquid



- ✓ Featuring a valve in the lid which prevents leakage if the mug is tipped over. Users can even drink (sucking) lying down
- ✓ Lids available with small aperture (most common), large aperture (for thickened liquids) and heat sensitive lid (changes colour if liquid too hot)
- ✓ Easy to hold - even with weak hands
- ✓ Suitable for warm and cold liquids
- ✓ Optional straw holder attachment holds the straw in place on the mug rim
- ✓ Choice of 4 colours
- ✓ Dishwasher, microwave and freezer safe
- ✓ CE marked under the Medical Device Directive 93/42/EEC



The drinking lid is a perfect fit for the Sure Grip Mug



The straw clip allows secure attachment of the drinking straw (not supplied) to the mug.



Thermo Lid, small aperture, indicates safe drinking temperature changing colour at 37° - 38°C



Specially designed drinking lids come with small or large aperture

Choose the mug and lid from the list below

OS FSA 5201	Suregrip Mug	Clear	200ml
OS FSA 5202	Suregrip Mug	Blue	200ml
OS FSA 5203	Suregrip Mug	Yellow	200ml
OS FSA 5204	Suregrip Mug	Red	200ml
OS FSA 5205	Suregrip Lid - Small Hole and Valve (most popular choice)		
OS FSA 5206	Suregrip Lid - Large Hole and Valve		
OS FSA 5207	Suregrip Thermo Lid - Small Hole and Valve		
OS FSA 5208	Suregrip Straw Clips		

Pack of 2

# PROVALE™ Cup

Uniquely, the Provale Cup delivers a set amount of fluid each time the cup is tipped. This important feature ensures a controlled delivery for those with swallowing difficulties.

- ✓ Designed with the assistance of interviews and feedback from over 500 speech-language therapists and occupational therapists
- ✓ Accurate, fixed delivery amount of 5cc and 10cc (depending on style) with each drinking motion
- ✓ Available with solid or clear base, blue lid (5cc) and brown lid (10cc)
- ✓ The clear base allows carer to see how much fluid has been delivered.
- ✓ Patented mechanism prevents over-delivery of fluid.
- ✓ Good nose clearance minimises head tilt
- ✓ Suitable for warm and cold liquids
- ✓ Can accommodate some thickened liquids
- ✓ Optional one or two easy-grip-handle or no-handle design
- ✓ Travel-mug appearance is functional, yet safe
- ✓ Easy cleaning microwave and dishwasher safe. No difficult valves or tubes to clean
- ✓ CE marked under the Medical Device Directive 93/42/EEC



Single handle clear



Accurate fixed delivery

Double handle clear

OS PDC 0501	PROVALE CUP	5cc	BLUE
OS PDC 0502	PROVALE CUP	5cc	BLUE/CLEAR
OS PDC 1001	PROVALE CUP	10cc	BROWN
OS PDC 1002	PROVALE CUP	10cc	BROWN/CLEAR


# NOVO Cup

The user is in complete control of the flow rate of fluid when the cap plug is in place by sucking on the spout. This cup is ideal for those with restricted movements of the head and neck.

- ✓ Suitable for warm and cold liquids
- ✓ Complete control of flow rate when cap plug is in place
- ✓ Ideal for those with restricted head and neck movement
- ✓ Spill proof when tipped over (with cap plug in place)
- ✓ Clear plastic beaker so contents are visible
- ✓ With measurement markings; contains 250ml liquid
- ✓ Saves time for users - no need to lift or support patients head whilst the patient is drinking
- ✓ Complete with two straws
- ✓ Dishwasher and microwave safe
- ✓ CE marked under the Medical Device Directive 93/42/EEC



OS FSA 5101	Novo Cup with Lid and 2 Straws	250ml
OS FSA 5102	Straws for Novo Cup	Pack of 50



## MAROON Spoons

Adults and children with feeding problems such as poor lip closure, tongue, throat and oral hypersensitivity, will find it easier to use the Maroon Spoons.

These sturdy spoons have shallow, narrow bowls. The design of the spoons makes it easier to get foodstuffs off the spoon - compared with traditional spoons.

- ✓ No sharp edges and produced from polycarbonate
- ✓ For hygiene reasons we recommend that Maroon Spoons be replaced regularly
- ✓ Shatter proof (not suitable for severe bite reflex)
- ✓ Lightweight
- ✓ Easy to clean and dishwasher safe
- ✓ Available in two sizes in packs of 20
- ✓ CE marked under the Medical Device Directive 93/42/EEC

**Dimensions:**

**Small Bowl**

Bowl dimensions, widest point: 25mm  
Bowl length: 35mm  
Length of spoon: 145mm

**Large Bowl**


Bowl dimensions, widest point: 31mm  
Bowl length: 40mm  
Length of spoon: 150mm

OS FSA 4901	Maroon Spoons Small	Pack of 20
OS FSA 4902	Maroon Spoons Large	Pack of 20

## DYSPHAGIA Cup

The Dysphagia Cup has been carefully designed, and is shaped to allow sufficient nose clearance - to empty contents without having to tilt the head back.

- ✓ Sized for efficient lip seal and control
- ✓ Shaped to direct flow of contents to the centre of the mouth
- ✓ Contoured for easier delivery of thickened liquid flow
- ✓ Suitable for warm and cold liquids
- ✓ Weighted to help reduce tremor frequency
- ✓ Choice of 2 colours, Almond or Clear Green. Clear cup allows carers to see how much liquid has been delivered
- ✓ One piece construction - more hygienic
- ✓ Durable, non toxic polycarbonate material
- ✓ Wide handle to allow extended palm inside so gripping not required
- ✓ Easy to clean, dishwasher safe
- ✓ CE marked under the Medical Device Directive 93/42/EEC




A successful swallow is one that has diverted food or liquid into the food tube (oesophagus) and not into the lungs. Frequently a person who has had a stroke, head injury or neurological disability has difficulty performing a successful swallow. Safe swallowing may be accomplished

by maintaining the chin in a tucked down position or by drinking thickened liquids or both.

The Dysphagia Cup has been designed to meet the needs of people who have difficulty in completing a successful swallow.

OS FSA 5002	Dysphagia Cup	200ml	Almond
OS FSA 5003	Dysphagia Cup	200ml	Clear Green



Kapitex have available Oral Motor Exercise tools for the rehabilitation of swallowing difficulties, particularly after stroke.




### ORALIGHT<sup>®</sup> Oral Motor Exercise System

Please contact Kapitex for a brochure and further information.



### DEPEND ON KAPITEX QUALITY AND SERVICE

Kapitex Healthcare provide a wide range of products dedicated to laryngectomy, tracheostomy, dysphagia and airway management. You can depend on our experience, our product quality - and our quality of aftersales service. To order any of the products in this brochure or for details of other Kapitex products contact Kapitex on 01937 580211, or by fax on 01937 580796 or write to us at the address below.



**Kapitex Healthcare Ltd.**  
Kapitex House, 1 Sandbeck Way, Wetherby,  
West Yorkshire LS22 7GH, England.  
Tel: 01937 580211 Fax: 01937 580796  
If dialling from outside the UK -  
Tel: +44 (0)1937 580211 Fax: +44 (0)1937 580796

© Kapitex Healthcare Ltd. 2004  
Due to continuous product developments it may be necessary to alter without notice products within our range. We endeavour to keep all customers informed of product improvements.  
P80S0011 03/04



ISO 9001 /2000  
EN 46001  
ISO 13485  
Cert. 257